



The MFD Automatic Polisher For Superior Polishing of Tablet Punches

With it

Reduce health & safety risk
Full GMP compliance
Operator time = 1 minute
Polish 60 punches in 15 minutes
Eliminate over-polishing
Polish the entire punch
Consistent results every time
Low skill level required
Minimum training needed
Improved finish reduces sticking/corrosion
Tools run smoother when highly polished

Typical return on investment 6 – 12 months

Web: www.casburt.com

Without it

Increase health & safety risk
Fails to meet GMP compliance
Operator time = several hours
Polish 60 punches in 180 minutes
Risk of over-polishing & damage
Polish the tip ONLY
Inconsistent results all the time
High skill level essential
High level training required
Poor finish increases sticking etc
Poorly maintained tools damage press

Ongoing high associated cost

Email: info@casburt.com

Suppliers of :- Tablet punches & dies – Tablet Presses – Spare Parts – Servicing - GMP Training